

## Product Information

### Solus Half Fraser Broth (FBH001, FBH002, FBH005)

<u>Formulation</u>	<u>Grams/litre</u>
Peptone mixture	15.0
Yeast extract	5.0
Aesculin	1.0
Disodium hydrogen phosphate	9.6
Potassium dihydrogen phosphate	1.35
Sodium chloride	20.0
Lithium chloride	3.0
Acriflavine	0.0125
Nalidixic Acid	0.01

### Directions

Weigh 55 grams of powder and dispense into 1 litre deionised water.

Allow to soak for 10 minutes, swirl to mix and sterilise at 121°C for 15 minutes.

Cool to 47°C and aseptically dispense into sterile containers.

Optional (not required for use with Solus Listeria test kits) : for visualisation of the aesculin reaction, prior to dispensing into final containers add 2 vials of F.A.C. supplement (SUP004).

### Appearance:

Straw opalescent broth with yellow fluorescence.



Certificate Number 11117  
ISO 9001

#### **Solus Scientific Solutions Ltd**

9 Mansfield Networkcentre, Millennium Business Park, Concorde Way, Mansfield, Nottinghamshire NG19 7JZ

**Tel** +44 (0)1623 429701 **Fax** +44 (0)1623 620977 **Email** enquiries@solusscientific.com

**Web** www.solusscientific.com

**Registered Address:** East Kilbride Laboratory, Torus Building, Rankine Avenue, East Kilbride, Glasgow, Lanarkshire G75 0QF Company Reg. No. SC361616, Registered in Scotland VAT No 982 2304 24