

Product Information

Solus Half Fraser Broth (FBH002, FBH005)

<u>Formulation</u>	<u>Grams/litre</u>
Peptone mixture	15.0
Yeast extract	5.0
Aesculin	1.0
Disodium hydrogen phosphate	9.6
Potassium dihydrogen phosphate	1.35
Sodium chloride	20.0
Lithium chloride	3.0
Acriflavine	0.0125
Nalidixic Acid	0.01

Directions

Weigh 55 grams of powder and dispense into 1 litre deionised water.

Allow to soak for 10 minutes, swirl to mix and sterilise at 121°C for 15 minutes.

Cool to 47°C and aseptically dispense into sterile containers.

Optional (not required for use with Solus Listeria test kits): for visualisation of the aesculin reaction, prior to dispensing into final containers add 2 vials of F.A.C. supplement (SUP004).

Appearance:

Straw opalescent broth with yellow fluorescence.

Solus Scientific Solutions Ltd

9 Mansfield Networkcentre, Millennium Business Park, Concorde Way, Mansfield, Nottinghamshire NG19 7JZ

Tel +44 (0)1623 429701 **Fax** +44 (0)1623 620977 **Email** enquiries@solusscientific.com

Web www.solusscientific.com

Registered Address: East Kilbride Laboratory, Torus Building, Rankine Avenue, East Kilbride, Glasgow, Lanarkshire G75 0QF Company Reg. No. SC361616, Registered in Scotland VAT No 982 2304 24