

# Certificate of Analysis



**Product Code:** EC-048S  
**Product:** Solus *E. coli* O157 ELISA (1 plate)  
**LOT / Batch:** EE032  
**Date of Manufacture:** 26.10.2020  
**Expiry Date:** 25.10.2021

This batch has been quality control tested to the following criteria:

Batch homogeneity CV%: < 10%.  
A panel of 8 non - *E. coli* O157 *Enterobacteriaceae* spp: Satisfactory  
Non-toxicogenic *E. coli* O157 NCTC 12900: Satisfactory

## Performance of In house Reference Standards

Sample	Net Signal OD	Specification OD
F2 <i>E. coli</i> O157 ~10 <sup>4</sup> cfu/ml	0.644	>0.300
F2 <i>E. coli</i> O157 ~10 <sup>5</sup> cfu/ml	1.503	>0.500
F2 <i>E. coli</i> O157 ~10 <sup>6</sup> cfu/ml	1.760	>1.000

Net Optical Densities read at 450nm of samples grown in mTSB + Novobiocin (20mg/L) for 20 hours at 41.5°C

	1	2
A	Negative control 0.000	<i>E5 Serratia liquifaciens</i> 0.044
B	Positive control 1.848	<i>B7 Salmonella Typhimurium</i> 0.016
C	Substrate blank 0.045	<i>A8 E. coli</i> 0.001
D	<i>F2 E. coli</i> O157 NCTC 12900 >2.500	<i>E2 E. coli</i> (Wild Type) 0.048
E	<i>A9 Proteus mirabilis</i> 0.044	Absorbances read at 450nm
F	<i>I3 Citrobacter freundii</i> 0.041	
G	<i>A10 Klebsiella pneumoniae</i> 0.028	
H	<i>E4 Shigella sonnei</i> 0.013	

# Certificate of Analysis



Production Batch is approved for release as it gives satisfactory QC performance.

Signed: Technical director/~~Finance Leader~~/Quality Manager

*D Schijwants*

Date:

3 | 11 | 20

