

# Certificate of Analysis

**Product Code:** SALSUPP-112.5  
**Product:** Solus One Salmonella Supplement  
**LOT / Batch:** SK10S  
**Date of Manufacture:** 13/12/2018  
**Expiry Date:** 31/12/2020

Satisfactory growth observed from the test organism in Buffered Peptone Water incubated for 24h at 37°C.

The supplement when incorporated into Mueller Hinton agar gave the expected growth/inhibition profile with 36 bacteria when incubated at 37°C for 24h.

Sterility of the supplement was confirmed as no growth was observed when added to blood agar incubated for 3 days at 37°C and 20°C.

<b><u>Performance:</u></b>	<b><u>Specification</u></b>	<b><u>Results</u></b>
<i>Salmonella typhimurium</i> NCTC 12023	Growth	Growth
<i>Staphylococcus aureus</i> ATCC 25923	No Growth	No Growth

The performance of this product is based on use of the supplement as detailed in the accompanying instructions for use.

**Gillian Bradley, Quality Manager**

**Date: 13/12/2018**

**Solus Scientific Solutions Ltd**

9 Mansfield Networkcentre, Millennium Business Park, Concorde Way, Mansfield, Nottinghamshire NG19 7JZ  
**Tel** +44 (0)1623 429701 **Fax** +44 (0)1623 620977 **Email** enquiries@solusscientific.com  
**Web** www.solusscientific.com

**Registered Address:** East Kilbride Laboratory, Torus Building, Rankine Avenue, East Kilbride, Glasgow, Lanarkshire G75 0QF Company Reg. No. SC361616, Registered in Scotland VAT No 982 2304 24