



Product Information

Solus Modified Peptone Water (MED038)

<u>Formulation</u>	<u>Grams/litre</u>
Tryptone	10.0
Sodium Chloride	5.0
Disodium Hydrogen Phosphate	3.6
Potassium Dihydrogen Phosphate	1.5
Neutraliser Mix	1.11

Directions

Weigh 21.21 grams of powder and disperse into 1 litre deionised water.
Allow to soak for 15 minutes, and swirl to mix. Autoclave at 121°C for 15 minutes.

Appearance:

Colourless to pale yellow

Solus Scientific Solutions Ltd

9 Mansfield Networkcentre, Millennium Business Park, Concorde Way, Mansfield, Nottinghamshire NG19 7JZ

Tel +44 (0)1623 429701 **Fax** +44 (0)1623 620977 **Email** enquiries@solusscientific.com

Web www.solusscientific.com

Registered Address: East Kilbride Laboratory, Torus Building, Rankine Avenue, East Kilbride, Glasgow,
Lanarkshire G75 0QF Company Reg. No. SC361616, Registered in Scotland VAT No 982 2304 24