

# Certificate of Analysis



**Product Code:** LIS1-0096  
**Product:** Solus ONE *Listeria* (1 plate)  
**LOT / Batch:** NDL019  
**Date of Manufacture:** 26.10.2020  
**Expiry Date:** 11.10.2021

This batch has been quality control tested to the following criteria:

Batch homogeneity: CV% < 10%.  
A panel of 9 non- *Listeria* spp: Satisfactory  
A panel of 11 *Listeria* spp: Satisfactory

## Performance of In house Reference Standards

Sample	Net Signal OD	OD Specification
<i>D8 L. monocytogenes</i> 4a LMG16783 5x10 <sup>4</sup>	>2.500	>0.500
<i>D8 L. monocytogenes</i> 4a LMG16783 1x10 <sup>5</sup>	>2.500	>1.500
<i>D8 L. monocytogenes</i> 4a LMG16783 1x10 <sup>6</sup>	>2.500	>2.000
<i>G3 L. innocua</i> wild type 5x10 <sup>4</sup>	2.001	>0.500
<i>G3 L. innocua</i> wild type 1x10 <sup>5</sup>	>2.500	>1.500
<i>G3 L. innocua</i> wild type 1x10 <sup>6</sup>	>2.500	>2.000
<i>D6 L. welshimeri</i> LMG11389 5x10 <sup>4</sup>	>2.500	>0.500
<i>D6 L. welshimeri</i> LMG11389 1x10 <sup>5</sup>	>2.500	>1.400
<i>D6 L. welshimeri</i> LMG11389 1x10 <sup>6</sup>	>2.500	>2.000

# Certificate of Analysis



Net Optical Densities of samples grown in SOLO+ at 30.0°C for 24 hours and boiled at 85-100°C for 20 minutes. Listeria cultures are diluted 1:100 in SOLO+ prior to analysis.

	1	2	3
A	Negative control 0.017	<i>D7 L. monocytogenes</i> >2.500	<i>E8 Brochothrix thermosphacta</i> 0.020
B	Positive control >2.500	<i>D8 L. monocytogenes</i> >2.500	<i>E9 Bacillus cereus</i> 0.022
C	Substrate blank 0.040	<i>D9 L. monocytogenes</i> >2.500	<i>E6 Enterococcus faecalis</i> 0.021
D	<i>G10 L. innocua</i> >2.500	<i>C10 L. monocytogenes</i> >2.500	<i>A8 Escherichia coli</i> 0.045
E	<i>D5 L. ivanovii</i> >2.500	<i>AH10 L. ivanovii londoniensis</i> >2.500	<i>C1 Staphylococcus aureus</i> 0.017
F	<i>D6 L. welshimeri</i> >2.500	<i>G3 L. innocua</i> >2.500	<i>E2 Escherichia coli (Wild Type)</i> 0.059
G	<i>D2 L. grayii</i> >2.500	<i>F1 Streptococcus bovis</i> 0.014	<i>F10 Enterococcus faecalis (Wild type)</i> 0.037
H	<i>D1 L. monocytogenes</i> >2.500	<i>E7 Kurthia zopfii</i> 0.041	Absorbances read at 450nm

Production Batch is approved for release as it gives satisfactory QC performance.

Signed: Technical director/Finance Leader/Quality Manager

Date:

9/11/20.