

Certificate of Analysis



Product Code: LIS1-0480
Product: Solus ONE *Listeria* (5 plate)
LOT / Batch: NDL026
Date of Manufacture: 21.06.2021
Expiry Date: 20.12.2022

This batch has been quality control tested to the following criteria:

Batch homogeneity: CV% < 10%.
A panel of 9 non- *Listeria* spp: Satisfactory
A panel of 11 *Listeria* spp: Satisfactory

Performance of In house Reference Standards

| Sample | Net Signal OD | OD Specification |
|----------------------------------------------------------|---------------|------------------|
| <i>D8 L. monocytogenes</i> 4a LMG16783 5x10 ⁴ | >2.500 | >0.500 |
| <i>D8 L. monocytogenes</i> 4a LMG16783 1x10 ⁵ | >2.500 | >1.500 |
| <i>D8 L. monocytogenes</i> 4a LMG16783 1x10 ⁶ | >2.500 | >2.000 |
| <i>G3 L. innocua</i> wild type 5x10 ⁴ | 1.046 | >0.500 |
| <i>G3 L. innocua</i> wild type 1x10 ⁵ | 1.595 | >1.500 |
| <i>G3 L. innocua</i> wild type 1x10 ⁶ | >2.500 | >2.000 |
| <i>D6 L. welshimeri</i> LMG11389 5x10 ⁴ | >2.500 | >0.500 |
| <i>D6 L. welshimeri</i> LMG11389 1x10 ⁵ | >2.500 | >1.400 |
| <i>D6 L. welshimeri</i> LMG11389 1x10 ⁶ | >2.500 | >2.000 |

Certificate of Analysis



Net Optical Densities of samples grown in SOLO+ at 30.0°C for 24 hours and boiled at 85-100°C for 20 minutes. Listeria cultures are diluted 1:100 in SOLO+ prior to analysis.

| | 1 | 2 | 3 |
|---|--------------------------------------|------------------------------------------------|-------------------------------------------------------|
| A | Negative control 0.027 | D7 <i>L. monocytogenes</i> >2.500 | E8 <i>Brochothrix thermosphacta</i> 0.007 |
| B | Positive control >2.500 | D8 <i>L. monocytogenes</i> >2.500 | E9 <i>Bacillus cereus</i> 0.045 |
| C | Substrate blank 0.043 | D9 <i>L. monocytogenes</i> >2.500 | E6 <i>Enterococcus faecalis</i> 0.003 |
| D | G10 <i>L. innocua</i> >2.500 | C10 <i>L. monocytogenes</i> >2.500 | A8 <i>Escherichia coli</i> 0.003 |
| E | D5 <i>L. ivanovii</i> >2.500 | AH10 <i>L. ivanovii londoniensis</i> >2.500 | C1 <i>Staphylococcus aureus</i> 0.023 |
| F | D6 <i>L. welshimeri</i> >2.500 | G3 <i>L. innocua</i> >2.500 | E2 <i>Escherichia coli</i> (Wild Type) 0.000 |
| G | D2 <i>L. grayii</i> >2.500 | F1 <i>Streptococcus bovis</i> 0.002 | F10 <i>Enterococcus faecalis</i> (Wild type) 0.005 |
| H | D1 <i>L. monocytogenes</i> >2.500 | E7 <i>Kurthia zopfii</i> 0.002 | Absorbances read at 450nm |

Production Batch is approved for release as it gives satisfactory QC performance.

Signed: Technical director/Operations Manager/Quality Manager

Date: 29/6/21