



Certificate of Analysis

Buffered Peptone Water pH (ISO) 7.0

Lot No: 202531

Product Code: MED017-10KG

Exp: 24/05/2022

Physical Characteristics

Dehydrated Medium: Fine Free Flowing Homogeneous/ Buff

Reconstituted Medium: Pale Straw

Precipitation: None

pH at 25°C: 7.1 pH specification: 7.0 +/- 0.2

Microbiological Performance

Organism Name	Organism Specification	Result Value
Salmonella typhimurium WDCM 00031	Recovery from 10-100cfu	Recovery from 10-100cfu
Salmonella enteritidis WDCM 00030	Recovery from 10-100cfu	Recovery from 10-100cfu
Escherichia coli WDCM 00013	Recovery from 10-100cfu	Recovery from 10-100cfu
Escherichia coli WDCM 00012	Recovery from 10-100cfu	Recovery from 10-100cfu
Escherichia coli WDCM 00013	+/- 30% colonies/T0	+/- 30% colonies/T0
Escherichia coli WDCM 00012	+/- 30% colonies/T0	+/- 30% colonies/T0
Staphylococcus aureus WDCM 00034	+/- 30% colonies/T0	+/- 30% colonies/T0
Listeria monocytogenes WDCM 00021	+/- 30% colonies/T0	+/- 30% colonies/T0
Listeria monocytogenes WDCM 00109	+/- 30% colonies/T0	+/- 30% colonies/T0

Product tested in compliance to ISO 11133:2014

Product is prepared for testing according to documented instructions. Cultures used for challenge testing are from approved culture collections and procedures used are industry accepted methodologies.

The information given above is believed to be correct, however, performance is only warranted when the product is prepared and used as described in the current version of the product specification.

Gillian Bradley – Quality Manager

Date: 24/05/2018

Solus Scientific Solutions Ltd

9 Mansfield Networkcentre, Millennium Business Park, Concorde Way, Mansfield, Nottinghamshire NG19 7JZ

Tel +44 (0)1623 429701 Fax +44 (0)1623 620977 Email enquiries@solusscientific.com

Web www.solusscientific.com

Registered Address: East Kilbride Laboratory, Torus Building, Rankine Avenue, East Kilbride, Glasgow, Lanarkshire G75 0QF Company Reg. No. SC361616, Registered in Scotland VAT No 982 2304 24