

Certificate of Analysis



Product Code: SOL001
Product: SOLO+ 500g
LOT / Batch: 212328
Date of Manufacture: 05/05/2020
Expiry Date: 04/05/2023

Physical Characteristics:

Dehydrated Medium: Fine free flowing homogenous/Buff
Reconstituted Medium: Dark straw with yellow fluorescence
Precipitation: None
pH at 25°C: 7.2
pH Specification 7.2 +/- 0.1

Microbiological Performance:

Organism Name	Organism Specification	Results Value
Listeria monocytogenes ATCC 19111	Recovery from 10-100cfu	Recovery from 10-100cfc
Listeria innocua NCTC 11288	Recovery from 10-100cfu	Recovery from 10-100cfc
Listeria ivanovii NCTC 11846	Recovery from 10-100cfu	Recovery from 10-100cfc
Listeria seeligeri NCTC 11856	Recovery from 10-100cfu	Recovery from 10-100cfc
Listeria welshimeri NCTC 11857	Recovery from 10-100cfu	Recovery from 10-100cfc
Enterococcus faecalis ATCC 29212	Recovery from 10-100cfc	Recovery from 10-100cfc
ELISA test performance	Satisfactory	Satisfactory

Product tested using methodologies in ISO 11133:2014

Product is prepared for testing according to documented instructions. Cultures used for challenge testing are from approved culture collections and procedures used are industry accepted methodologies.

The information given above is believed to be correct, however, performance is only warranted when the medium is prepared and used as described in the instructions accompanying the goods either in or on the packaging.

Gillian Bradley, Quality Manager

Date: 28/05/2020