

Certificate of Analysis



| | |
|-----------------------------|-----------------------------|
| Product Code: | SUP011 |
| Product: | P.A.C. Selective Supplement |
| LOT / Batch: | 457876 |
| Date of Manufacture: | 31/01/2021 |
| Expiry Date: | 31/01/2023 |

Quality Assurance

Products are manufactured in accordance with an ISO9001/13485 Quality system which covers design, manufacture, and supply of in-vitro diagnostic devices and consumables.

Quality Control Testing

The above products were tested in the Quality Control Laboratory and met defined acceptance criteria using our standard test method. The Quality Control Laboratory method complies with BS EN 12322:1999 "In vitro diagnostic medical devices - Culture media for microbiology - Performance criteria for culture media".

Physical Characteristics: Appearance of both lyophilised and reconstituted product was satisfactory.

Performance: Tests on content, microbiological performance with a range of challenge organisms and bioburden were within specification.

Gillian Bradley, Quality Manager

Date: 31/01/2021